
# SOCOREX® classic 167/187/187F/287

## ANIMAL HEALTH SYRINGES VETERINÄR SPRITZEN SERINGUES VÉTÉRINAIRES

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#### ⚠ SAFETY PRECAUTIONS

- Read the operating instructions carefully before use. Keep available for future reference
- Observe the manufacturer's recommendations
- Follow instructions for proper use of injectables
- Before each use, check instrument for perfect working condition
- Purge, rinse and dry the instrument carefully after each use
- Never try to separate blocked parts by using force
- Repeat use of the same needle may carry over contaminations to other animals
- Instrument manufacturer cannot be made responsible for needle breakage during use
- Mind possible hand fatigue and its medical consequences
- Do not use syringe for purposes other than intended for
- A saline environment may hasten replacement of instruments parts
- Use only original spare parts

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### DESCRIPTION

Congratulations for buying a fine quality Swiss product. The self-refilling syringes 167, 187, 187F and 287 with spring loaded plunger and three-way valve system are intended for mass vaccination of poultry, swine, cattle, fish and small laboratory animals. All models can be autoclaved fully assembled at 121°C (20 minutes). These automatic syringes allow trouble free injection of volumes as small as 0.025 mL, various water or oil based vaccines or difficult heavy iron solutions with an excellent reproducibility. Valves are fitted with stainless steel balls for improved reliability. Volume is easily set with a micrometric screw and a lock nut. The graduated ground glass barrel is interchangeable.

- Model **167 with two ring handle**, ranging from 0.025 to 5 mL. Supplied with feed tubing but compatible with vial holder.
- Model **187 with pistol handle**, ranging from 0.025 to 10 mL. Supplied with feed tubing, vial holder or both.
- Model **187F with pistol handle**, with volumes 0.2, 0.3, 0.5 and 1 mL. Supplied either with feed tubing or vial holder. Other volumes available on request
- Twin model **287 with pistol handle**, ranging from 0.025 to 5 mL. Supplied with feed tubing and with vial holder, they are efficient productivity booster by allowing injection of two different compounds at the same time.

Veterinary needles and drenchers are to be ordered separately. Ask for more details about our large selection of Luer Lock stainless steel needles.

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## HOLDING THE SYRINGE

Ergonomic, user friendly handle, perfect balance and smooth plunger travel assure fatigue-free operation even during mass vaccinations at rates of up to 800 injections per hour.

To operate 167 model, place forefinger and middlefinger in the two rings of the syringe handle and move spring-loaded plunger with the thumb.

On 187, 187F and 287 models, place hand on handle and move plunger with the palm of the hand.

To change position of the syringe handle (41), untighten screw (44), turn syringe body to the desired position and retighten screw (44).

## OPERATION

All self-refilling syringes are supplied ready to use, no assembling or preparation is required.

### • Self-refilling syringe with vial holder

To draw liquid out of an original vaccine bottle with septum closure, remove lock (55), attach vaccine vial (penicillin neck, Ø 20 mm) to intake needle (52), slide lock in again. Rubber washer (1.187.745) improves tightness. Disassemble vial holder by removing lock (55) and ring (54).

### • Self-refilling syringe with feed tubing

To aspirate liquid from a bottle with rubber stopper carried at the belt, connect one end of tube (23) to the syringe. Puncture bottle cap with feed needle (26) and vent cannula (25). Connect feed tube to the feed needle. When the bottle is hanging, with neck towards the bottom, it is recommended to use the optional draw off needle with air inlet and filter (ref. 1.187.590).

To aspirate liquid from an open vial, connect sinker (24) to the end of feed tube (23) and immerse in liquid.

PP reservoir of 60, 125 or 250 mL capacity (see list of accessories) can be fitted on all syringe types, after adaptation of a Rotlauf intake valve (ref. 187.613).

Attach injection needle (not supplied) to nozzle (14). Optional manifold (ref. 1.287.685) allows injection through one needle on Twin 287 model.

## Volume setting (except 187F)

Loosen locking nut (4) and adjust plunger edge to the desired volume graduation on glass barrel by turning socket (3). Tighten locking nut (4) in order to prevent any volume alteration during plunger activation.

On 287 Twin syringe, the volume is independently adjustable on each barrel.

## Priming

Prime valve system by activating plunger until liquid fills up the syringe. Set syringes of small capacity on largest volume to prime. To expel remaining air bubbles, hold the syringe upright and give the plunger a few short strokes. The liquid lost during priming may be collected in a clean recipient and returned to the reservoir.

## End of dosing

When dosing is completed, remove feed sources and activate plunger several times in order to eliminate any residual liquid. Rinse, disassemble, clean and dry the instrument parts according to the instructions in the "maintenance" chapter. **Never store a syringe without prior purging, cleaning and drying.**

## Manifold mounting (Twin model 287)

Loosen both locking nuts (4) and both handle screws (44). Position the manifold on both Luer lock nozzles (14). Tighten both manifold security rings. Tighten both locking nuts, then both handle screws.

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## MAINTENANCE

### Cleaning

The self-refilling syringe and accessories must be carefully cleaned and dried after each use.

Therefore, aspirate an appropriate cleaning solution (water or detergent) and rinse thoroughly. Activate the syringe several times with air to remove residual liquid before disassembling.

### Disassembling

Remove needle, tubing and/or vial holder. Hold handle (41) and unscrew casing (8) in order to access the glass barrel (7). Push rod (1) or handpiece (48) (Twin 287) to free plunger (6). Unscrew aspiration (50) and dispensing valves (14). Dry the valve assembly (seat, ball), syringe body, barrel and plunger before reassembling and storing.

Using on a regular basis the ABF Premium lubricant after cleaning/drying will protect the plunger and extend the syringe lifespan.

### Trouble shooting

A clogged aspiration valve can easily be fixed up by gently detaching the ball from the valve seat. Therefore, introduce a small pin into the intake valve opening and push until the ball moves. If this operation is not possible, immerse the valve in hot water containing an appropriate detergent and operate the plunger. Make sure the parts are rinsed thoroughly and dried completely.

To prevent a mispositioning of aspiration valve ball, do not release spring loaded plunger rod abruptly while aspirating.

During the distribution, a thin film of liquid between the plunger and the barrel guarantees smooth plunger movement and tightness. After long series of injections, some liquid may be found in the back of the plunger and should be removed when vaccination or dosing is completed.

Only use original spare parts to replace any damaged components of the syringe.

### Sterilization

*Warning: optional PVC protection sleeve on glass barrel and vial holder rubber washer are not autoclavable. If used, remove them before autoclaving the instrument. Do not sterilize a syringe that has not been previously primed and cleaned. Slightly untighten the syringe body and valves before sterilization. The self-refilling syringe can be boiled in water at 100°C during 30 minutes. It is also possible to sterilize in an autoclave at 121°C during 20 minutes (1 Bar). Retighten all connections after sterilization and insure the instrument is dried during prolonged period of time. Correct sterilization conditions are the responsibility of the user.*

## MATERIALS

All parts in contact with the liquid are made of carefully selected materials:

- Feed tubing: silicone
- Barrel: sodocalcic AR glass
- Plunger: stainless steel
- Barrel washer: FPM
- Valve washers: PTFE
- Valve springs and balls: stainless steel
- Syringe casing, nozzle valve, intake valve, sinker and aspiration cannula: nickel/chrome plated brass

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## PERFORMANCE AND ORDERING INFORMATION

### Package contents

Standard supply includes Luer Lock nozzle, set of spares for valves (springs, valve balls and washers), operating instructions and according to the model:

- 187.1, 187.1.F and 287.1: vial holder set (vial holder, locking ring, bottle lock, tube connector, feed needle, rubber washer).
- 167.2, 187.2, 187.2.F and 287.2: feed tube set (silicone tubing, sinker, aspiration and vent cannulas). Twin 287 includes two each of the above accessories.
- 187.3: both vial holder and feed tube sets.
- 187: 2 mL and 5 mL, additional spring for viscous liquids.

Optional threaded Rotlauf and Record Unifix nozzles available on request.

Work station to be ordered separately.

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
<sup>1)</sup> at 0,075 mL

167.2 model with feed tube Cat. No.	187.1 model with vial holder Cat. No.	187.2 model with feed tube Cat. No.	187.3 model with vial holder and feed tube Cat. No.	287.1 Twin model with vial holder Cat. No.	287.2 Twin model with feed tube Cat. No.
167.2.05003	187.1.05003	187.2.05003		287.1.05003	287.2.05003
167.2.05005	187.1.050053	187.2.05005	187.3.05005	287.1.05005	287.2.05005
167.2.0501	187.1.0501	187.2.0501	187.3.0501	287.1.0501	287.2.0501
167.2.0502	187.1.0502	187.2.0502	187.3.0502	287.1.0502	287.2.0502
167.2.0505	187.1.0505	187.2.0505	187.3.0505	287.1.0505	287.2.0505
	187.1.0510	187.2.0510	187.3.0510		

Volume* (mL)	Innaccuracy (E%)	Imprecision (CV%)	187.1.F model with vial holder Cat. No.	187.2.F model with feed tube Cat. No.
0.2	< +/- 2.5%	< 0.5%	187.1.F05002	187.2.F05002
0.3	< +/- 2.0%	< 0.4%	187.1.F05003	187.2.F05003
0.5	< +/- 1.5%	< 0.3%	187.1.F05005	187.2.F05005
1	< +/- 1.5%	< 0.3%	187.1.F0501	187.2.F0501

\* Other volumes available upon request.

### Viscosity data

Variable Model mL	Maximal cinematic viscosity* (cst)	Example of liquids
0.3	1800	SAE 50 oil
0.5	2500	Fluid honey
1	1200	Silicon oil

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Variable Model mL	Maximal cinematic viscosity* (cst)	Example of liquids
2	400 (up to 600) <sup>1)</sup>	Glycerol
5	150 (up to 200) <sup>1)</sup>	SAE 20 oil
10	20	Cream

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Test conditions to determine max. viscosity
Syringe with horizontal feed tubing Syringe set at maximal volume Aspiration time < 3 sec.

Fixed Model mL	Maximal cinematic viscosity* (cst)	Example of liquids
0.2	2500	Fluid honey
0.3	2500	Fluid honey
0.5	2500	Fluid honey
1	1200	Silicon oil

Test conditions to determine max. viscosity
Syringe with horizontal feed tubing Aspiration time < 3 sec.

\* Cinematic viscosity in centistokes (cst). Dynamic viscosity in centipoises (cps) obtained by multiplying cinematic viscosity by liquid density (kg/dm<sup>3</sup>).

<sup>1)</sup> When using viscous liquid spring (see accessories)

## ACCESSORIES

Cat. No.	Description	Packaging
<b>PVC barrel protection sleeve</b>		
1.170.4005	0.2 / 0.3 / 0.5 mL	6 / pk
1.170.401	1 mL	6 / pk
1.170.402	2 mL	6 / pk
1.170.405	5 mL	6 / pk
1.170.410	10 mL	6 / pk
<b>PP reservoir with Rotlauf thread</b>		
1.187.060	60 mL	1 pc
1.187.150	125 mL	1 pc
1.187.250	250 mL	1 pc
<b>Drenchers – mouth &amp; nose cannula</b>		
376.0925S	Straight, 25 mm	3 / pk
376.1251S	Straight, 51 mm	3 / pk
376.16102S	Straight, 102 mm	2 / pk
376.30203S	Straight, 203 mm	2 / pk
376.0925B	Bent, 25 mm	3 / pk
376.1251B	Bent, 51 mm	3 / pk
376.16102B	Bent, 102 mm	2 / pk
376.30203B	Bent, 203 mm	2 / pk
<b>Optional one needle manifold</b>		
1.287.685	For 287 Twin syringe	1 pc
<b>Hard extension tubing</b>		
1.187.20	20 cm	1 pc
1.187.35	35 cm	1 pc

### Fish injections

FishGuide™ is a convenient accessory providing for correct needle positioning while injecting fishes. Fitting Socorex® syringes 167/187 models, it is most helpful in all fish farms. The accessory is supplied in three sizes: 15 mm, 25 mm and 45 mm and do not require any tool to be installed and removed. It can be positioned over 360° on the syringe. The two large copper grips facilitate fish positioning and their angle can be enlarged according to fish size.

**Notes:** Improper FishGuide™ positioning or overtightening screw may damage nozzle and impair needle fitting. Reduce risk of injuries by removing FishGuide™ when changing needle.

#### Assembling (Fig. 1, 2 and 3)

1. Attach needle (Fig. 1).
2. Insert FishGuide™ until abutment, adjust position, softly tighten thumb screw (Fig. 2).
3. Adjust angle as necessary (Fig. 3).


Fig. 1

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
Fig. 3

## WARRANTY

The self-refilling syringe is guaranteed for a period of one year against any manufacturing or material defect. Cases due to non-respect of manufacturer's instructions are excluded from the warranty. Glass barrel breakage or damages to the plunger due to improper handling or cleaning are not covered by the warranty.

Should regular maintenance not eliminate a detected defect, return instrument to the dealer from whom it was purchased, after having obtained a return authorisation.

**Rinse, dry and decontaminate instrument prior to its return.**

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