

# 02	187.751	Spring, rod		0.3 to 1 mL
	832.751			2 mL
	187.752	Viscous liquid spring, rod		2 mL
	832.752	Spring, rod		5 mL
	187.753	Viscous liquid spring, rod		5 mL
	187.754	Spring, rod		10 mL
	187.755			20 mL
# 03	187.961	Socket (incl. part 46)		0.3 to 1 mL
	187.962			2 mL
	187.963			5 mL
	187.964			10 mL
	187.965			20 mL
# 04	187.551	Front nut		0.3 to 1 mL
	187.552			2 mL
	187.553			5 / 10 mL
	187.554			20 mL
# 05	187.520	Cap		0.3 / 0.5 mL
	187.521			1 mL
	187.522			2 mL
	187.523			5 mL
	187.524			10 mL
	187.525			20 mL
# 06	187.584	Plunger, piston	Stainless steel	0.3 mL
	187.585			0.5 mL
	187.586			1 mL
	187.587			2 mL
	187.588			5 mL
	187.589			10 / 20 mL
# 07	0187.003	Barrel (6 pces)	Glass	0.3 mL
	0187.005			0.5 mL
	0187.01			1 mL
	0187.02			2 mL
	0187.05			5 mL
	0187.10			10 mL
	0187.20			20 mL
# 08	187.502	Casing	Chrome plated brass	0.3 / 0.5 mL
	187.503			1 mL
	187.505			2 mL
	187.507			5 mL
	187.509			10 mL
	187.510			20 mL
# 12	187.980	Spring valve (10 pces)	<i>Available as a set</i>	0.3 mL
	187.981			0.5 to 10 mL
	187.982	Spring, intake valve (10 pces)		20 mL
	187.983	Spring, outlet valve (10 pces)		20 mL
# 13	187.921	Washer, PTFE, (10 pces)	<i>Available as a set</i>	0.3 to 10 mL
	187.922			20 mL

# 14	705.010	Outlet nozzle valve, Luer Lock	Chrome plated brass	0.3 to 10 mL	
	707.010	Outlet nozzle valve, Record Unifix		0.3 to 10 mL	
	714.010	Outlet nozzle valve, Luer Linden w/o locking nut		0.3 to 10 mL	
	714.910	Outlet nozzle valve, Luer Linden with locking nut		0.3 to 10 mL	
	720.010	Outlet nozzle valve, Rotlauf		0.3 to 10 mL	
	186.562	Outlet nozzle valve for Needle Guard (Sterimatic®), Luer Lock		0.3 to 10 mL	
	705.020	Outlet nozzle valve, Luer Lock		20 mL	
# 15	187.928	Washer (10 pces)	<i>Available as a set</i>	FPM	0.3 / 0.5 mL
	187.924				1 mL
	187.925				2 mL
	187.926				5 mL
	187.927				10 / 20 mL
# 17	187.991	Valve ball (10 pces)	<i>Available as a set</i>		0.3 to 10 mL
	187.993				20 mL
# 23	187.705	Feed tubing	1 meter / pk		0.3 to 10 mL
	187.910		10 meter / pk		0.3 to 10 mL
	187.710		by the meter		0.3 to 10 mL
	187.707		1 meter / pk		20 mL
# 24	187.596	Sinker	Chrome plated brass		0.3 to 10 mL
# 25	187.933	Vent cannula (10 pces)	Plastic		0.3 to 5 mL
	187.591	Vent cannula	Metal		10 mL
	187.593				20 mL
# 26	187.592	Feed needle	Chrome plated brass		0.3 to 10 mL
	187.594				20 mL
# 41	187.560	Handle, blue			0.5 to 10 mL
	187.570				20 mL
# 44	187.701	Screw			0.5 to 10 mL
	187.932	Screw and brake			20 mL
# 46	187.711	Sleeve, PTFE			0.3 to 1 mL
	187.712				2 mL
	187.713				5 mL
	187.714				10 / 20 mL
# 50	187.959	Intake valve, ball type	Chrome plated brass		0.3 to 1 mL
	187.960		Wide hole size, serial Nos. > 13021001		2 to 10 mL
	187.616	Intake valve, special outlet for 5 mm tubing			0.3 to 10 mL
	187.613	Intake valve, Rotlauf threaded			0.3 to 10 mL
	187.621	Intake valve			20 mL
# 51	188.903	Connecting tube for vial holder model, silicone (10 pces)			0.3 to 10 mL
# 52	188.952	Intake needle (incl. part 51)			0.3 to 10 mL
# 53	187.950	Vial holder (2 faces), blue			0.3 to 10 mL
# 54	188.513	Ring, vial holder			0.3 to 10 mL
# 55	187.740	Lock, vial holder, blue			0.3 to 10 mL
# 56	1.187.745	Tightness rubber washer, Ø 20 mm (10 pcs)			0.3 to 10 mL
# 70	187.935	Washer, Rotlauf (10 pces)			0.3 to 10 mL
# 72	188.556	Vial holder with feed needle, Rotlauf			0.3 to 10 mL